



IGG (HIZENTRA) CONTINUOUS SUBCUTANEOUS INFUSION

PROCEDURE FOR SQ INFUSION:

- A. Wash hands
- B. Gather supplies:
 1. Chloraprep, swab sticks (3)
 2. Sub Q needle sets with transparent dressing included (type determined by order: i.e. single vs. multiple)
 3. Injection cap
 4. Tape
 5. Freedom 60 pump
 6. Freedom 60 tubing (rate as per order)
 7. 60cc syringe luer lock
 8. 20g needle x 1
 9. Sharps
 10. 10ml saline flushes
 11. 1 box of alcohol wipes
- C. Attach the injection cap to the end of the SQ set
- D. Attach the 5ml syringe of saline to the injection cap and prime the SQ set
- E. Choose appropriate sites for individual patient
 - a. Appropriate insertion sites include: abdomen, upper buttocks, lateral thigh or lateral hip
 - b. Inappropriate insertion sites are: at the waistline, near bony prominences and anywhere there is not intact skin
- F. Cleanse each site with a chloraprep swab using a back and forth friction scrub for at least 30 seconds. Allow area to dry. DO NOT BLOT.
- G. Insert subcutaneous needle(s) and secure with transparent dressing
- H. Check for a blood return at sites by gently pulling back on saline syringe. If there is a blood return, remove Sub Q needle and begin again. If no blood return noted, proceed to next step.



IGG CONTINUOUS SUBCUTANEOUS INFUSION

- I. Mix IGG per plan of treatment and draw up into the 60cc syringe
- J. Attach tubing labeled “Freedom 60” to 60cc syringe then place syringe into pump
- K. With syringe inserted into the Freedom 60 pump, prime the Freedom 60 tubing
- L. Wipe the injection cap with an alcohol wipe using a vigorous scrub for 30 seconds and let air dry for 60 seconds
- M. Attach end of Freedom 60 tubing to the injection cap
- N. Turn on the pump to begin infusion
- O. Secure tubing with tape
- P. When infusion is complete, turn the pump to the off position
- Q. Remove SQ needles and cover the site with a band aid
- R. Discard needle set, pump tubing and syringe into the sharps container
- S. Wash hands

IMPORTANT SAFETY ALERT:

TO AVOID RAPID INFUSION OF YOUR MEDICATION, OR AN ERROR, PLEASE VERIFY THE TUBING IS LABELED “FREEDOM 60” BEFORE INFUSING YOUR MEDICATION.

Revised 8/2015

